

Caspar™ Lumbar Retractor System



Aesculap Spine

Caspar™ Lumbar Retractor System



Product Description

The CASPAR Lumbar Retractor System (CLR) is ideal for medial and paraspinous approaches to perform discectomies, pedicle screw fixation, and many other lumbar procedures. The CLR is a stable system for safely retracting muscles and other tissue during surgery, allowing for optimal visualization.

Blade Widths and Lengths Available:
20 mm for monosegmental access
40 mm for bi-segmental access
60 mm for multi-segmental access



Fenestrated titanium blades provide excellent view during radiographic imaging.



Optimal tooth design provides excellent anchorage while preserving soft tissue.



Retractor hinge improves positioning and provides quick, easy and precise handling characteristics to ensure comfort for the surgeon



Side loading, quick locking mechanism makes the assembly of the retractors and blades quick and easy.

Caspar™ Lumbar Retractor System

Blade Options

20 mm Blades for Mono-segmental Access

Item No.	Description
BW020T	20 x 35 mm
BW021T	20 x 40 mm
BW022T	20 x 45 mm
BW023T	20 x 55 mm
BW024T	20 x 65 mm



40 mm Blades for Bi-segmental Access

Item No.	Description
BW030T	40 x 35 mm
BW031T	40 x 40 mm
BW032T	40 x 45 mm
BW033T	40 x 55 mm
BW034T	40 x 65 mm
BW035T	40 x 75 mm







60 mm Blades for Multi-segmental Access

Item No.	Description
BW040T	60 x 35 mm
BW041T	60 x 40 mm
BW042T	60 x 45 mm
BW043T	60 x 55 mm
BW044T	60 x 65 mm
BW045T	60 x 75 mm
BW046T	60 x 85 mm



Instrument Overview

Item No.	Description	
BW028T	Hemi-laminectomy Hook, Small	
BW001R	Lumbar Retractor	
BW011P	Depth Gauge Instrument	
BW038T	Hemi-laminectomy Hook, Large	
BW010R	Blade Handle	
<p><i>Note: Two handles are required for each system.</i></p>		
TE749R	Hex Key	
BW121P	CLR Tray	

Caspar™ Lumbar Retractor System

Recommended Set Configuration and Order Information

Item No.	Description	Recommended Qty.
Instruments		
BW001R	CLR System, Retractor ONLY	1
TE749R	Hex Key	1
BW010R	Blade Handle	2
BW011P	Depth Gauge Instrument	1
BW028T	Hemi-laminectomy Hook, Small	1
BW038T	Hemi-laminectomy Hook, Large	1
Blades		
BW020T	Retractor Blade, 20 x 35 mm	2
BW021T	Retractor Blade, 20 x 40 mm	2
BW022T	Retractor Blade, 20 x 45 mm	2
BW023T	Retractor Blade, 20 x 55 mm	2
BW024T	Retractor Blade, 20 x 65 mm	2
BW030T	Retractor Blade, 40 x 35 mm	2
BW031T	Retractor Blade, 40 x 40 mm	2
BW032T	Retractor Blade, 40 x 45 mm	2
BW033T	Retractor Blade, 40 x 55 mm	2
BW034T	Retractor Blade, 40 x 65 mm	2
BW035T	Retractor Blade, 40 x 75 mm	2
BW040T	Retractor Blade, 60 x 35 mm	2
BW041T	Retractor Blade, 60 x 40 mm	2
BW042T	Retractor Blade, 60 x 45 mm	2
BW043T	Retractor Blade, 60 x 55 mm	2
BW044T	Retractor Blade, 60 x 65 mm	2
BW045T	Retractor Blade, 60 x 75 mm	2
BW046T	Retractor Blade, 60 x 85 mm	2
Trays		
BW121P	CLR Tray	1
Recommended Sterile Containers		
JN441	Full-size Perf. Bottom, 5 1/3"	1
JK485	Full-size Lid, with Retention Plate (Red)	1
JG786U	Identification Label, Blue	2

Set No.	Set Description
ST0156	Caspar Lumbar Retractor (CLR) Set

For additional information, contact Customer Service at 1-866-229-3002 or visit our website at www.aesculapimplantsystems.com.

All rights reserved. Technical alterations are possible. This leaflet may be used for no other purposes than offering, buying and selling of our products. No part may be copied or reproduced in any form. In the case of misuse we retain the rights to recall our catalogs and price lists and to take legal actions.

©2016 AESCULAP. ALL RIGHTS RESERVED. PRINTED IN THE USA.
Aesculap is an equal opportunity employer

Aesculap Implant Systems, LLC | 3773 Corporate Parkway | Center Valley, PA | 18034
Phone 866-229-3002 | Fax 610-984-9096 | www.aesculapimplantsystems.com

Aesculap Implant Systems, LLC - a B. Braun company

DOC530 Rev B 11/16